



# POSEIDON

## SINONASAL IRRIGATION/SUCTION SYSTEM

*See with Clarity, Operate with Confidence - Simplified for Every Sinus Procedure*

### Ergonomic Surgeon-Centric Design: Available in Mini, Small, Medium, & Large

- Reduces fatigue and improves grip to maximize procedural comfort and control



### WaveGrip™ Sheath secures to endoscope, maintaining spatial orientation

- Provides quick confirmation of septal landmarks reducing procedural time

### Detachable Field Irrigation Feature for Targeted Fluid Delivery

- Quick-connect irrigation port allows immediate access for targeted tissue irrigation without interrupting the sterile field.

### Hands-Free Irrigation/Suction Control

- Pneumatic foot pedals activate separate fluid and vacuum channels for enhanced ease of control

### Valve Cartridge - Smooth Operation, Controlled Delivery

- Consistent Irrigation delivery and independent suction control
- Independent control of irrigation and suction

Product Code	Description	Qty/CS
SCI-LV-MINI	WaveGrip™ Sheath & Line Set –MINI	6
SCI-LV-S	WaveGrip™ Sheath & Line Set -SMALL	6
SCI-LV-M	WaveGrip™ Sheath & Line Set -MEDIUM	6
SCI-LV-L	WaveGrip™ Sheath & Line Set -LARGE	6
SCI-POTOW	Poseidon Control Tower	1
SCI-FTSWTCH	Poseidon Replacement Control Pedal	1
SCI-CUFF	Poseidon Replacement IV Pressure Cuff	1



Huck Medical Technologies, Inc.  
 111 Cash Street  
 Jacksonville, Tx 75766  
 903-586-4555  
 www.huckmedtech.com  
 customerservice@huckmedtech.com



LIT-SCI-PLC-010 REV 3.0.0